English

Fig 1 Shaker Incubator controls




Fig 2 Platforms and shaking actions

SI50: Orbital shaker
(includes cradle system with four securing bars)

SI70: Gyrotory rocker
(includes blue non-slip mat and tier system)

SI80: See-saw rocker
(includes blue non-slip mat and tier system)

Thank you for purchasing this Barloworld Scientific product. To get the best performance from the equipment, and for your personal safety, please read these instructions carefully before use.

This equipment is designed to operate under the following conditions:-

- ❖ For indoor use only
- ❖ Use in a well ventilated area
- ❖ Ambient temperature range +5°C to +40°C
- ❖ Altitude to 2000m
- ❖ Relative humidity not exceeding 80%
- ❖ Mains supply fluctuations not exceeding 10%
- ❖ Overvoltage category II IEC60364-4-443
- ❖ Pollution degree 2 IEC664

If the equipment is not used in the manner described in this manual the protection provided by the equipment may be impaired.

Electrical installation

THIS EQUIPMENT MUST BE EARTHED

Before connection please ensure that the line supply corresponds to that stated on the rating label.

Power requirements all models:

Voltage & Wattage
230V, 50Hz, ~ , 350W

There is an IEC socket at the rear of the instrument for connection to the mains supply. The unit is supplied with two mains leads fitted with IEC plugs for connection to the instrument. One has a U.K. 3 pin plug and the other has a 2 pin "Shuko" plug for connection to the mains. Choose the lead appropriate for your electrical installation and discard the other.

Should neither lead be suitable, take the lead with the U.K. plug and replace the plug with a suitable alternative. See the enclosed instruction sheet for advice on how to carry out this procedure.

Should the mains lead require replacement a cable of 1mm² of harmonised code H05W-F connected to an IEC 320 plug should be used.

IF IN DOUBT CONSULT A QUALIFIED ELECTRICIAN

The mains lead should be connected to the instrument BEFORE connection to the mains supply.

Safety Advice Before Use

- ❖ Never move or carry the unit when in use or connected to the mains electricity supply.
- ❖ Do not mix flammable liquids or use the equipment in hazardous atmospheres.
- ❖ In the case of mains interrupting, a fault or mechanical failure, the unit will continue to operate on removal of fault.
- ❖ Mechanical energy can lead to breakage of glass vessels. Use with care.

General description

CAUTION: Keep fingers clear when operating door.

This range of Stuart incubators contains three models each fitted with a moving platform offering a choice of shaking action.

SI50: Has an orbital shaker with variable speed between 20 and 250rpm and an orbit of 16mm. The shaker platform is fitted with an robust mat and adjustable rubber coated bars to secure the load.

SI70: Has a 3D gyratory rocker with variable speed between 10 and 70rpm and an angle of tilt adjustable from 3-12°. The shaker platform is fitted with an anti-slip mat.

SI80: Has a simple side to side rocking platform with variable speed between 10 and 70rpm and an amplitude adjustable from 1-7°. The shaker platform is fitted with an anti-slip mat.

The incubator compartment is heated by an efficient sheathed element. Forced air circulation and precise digital control gives accurate heat distribution and a quick response to set temperatures. For visibility of the contents the walls of the incubator are made from smoked acrylic.

The moving platform is operated by a geared induction motor. The speed is modulated by a soft start feedback controller in order to avoid splashing on start up and maintain long-term stability of set speed.

SI70 and SI80 are supplied with 2 extra platforms that may be mounted in a tier to increase capacity.

The doors of all models are fitted with a safety interlock device to prevent the platform moving if the door is open.

Operation (see fig. 1 & 2 for controls)

Position the apparatus on a firm level surface ensuring there is sufficient room to open the door outward.

Ensure the temperature and speed control knobs are fully to the left and switch the mains, motor and heater ON as required.

Temperature Selection - all models

Depress and hold the press to set switch and turn the temperature control knob until the required temperature appears on the display. On releasing the switch, the display will revert to the actual temperature of the incubator. The load indicator will remain lit until the desired temperature setting is reached at which time it will flash to show that the temperature is being controlled.

Once the temperature has been set, turn the over temperature control knob to approximately 20°C above the set temperature.

Platform operation - SI50

Secure the vessels to be shaken between the rubber coated bars. Ensure the door is fully closed and then turn the speed controller to the right until the required shaking speed is observed.

Platform operation - SI70 and SI80

Place the vessels to be shaken on the platform with the load distributed as evenly as possible.

Ensure the door is fully closed and then turn the speed controller to the right until the required shaking speed is observed.

Installing the tier system - SI70 and SI80

SI70 and SI80 are supplied with 2 extra platforms tiers which may be used to increase capacity. The units are supplied with 4 long magnetic support legs and 8 shorter magnetic support legs to allow a choice of spacing between the platforms. Also both platforms may be fitted at the same time if the shorter legs are used.

To fit the accessory platform, first screw a magnetic leg onto the captive thread in each corner.

Ensure that the motor is turned off. The accessory platform can then simply be placed on top of the main platform and will be held in place by the magnetic legs. The second accessory platform can then be similarly located on top of the first. To remove for cleaning etc., simply lift the platform, separating the magnetic legs from the supporting platform below.

Adjustment of Tilt Angle - SI70

1. Switch motor off.
2. Remove anti-slip mat from the platform and switch the motor back on and adjust the motion to a slow speed. Observe the two apertures at the centre of the platform, during rotation two adjusting screws will come into view. When they are aligned switch the motor off.
3. There are 2 locking screws located under the platform on the drive boss, facing front on. Check that these are accessible and, if not, switch the motor back on and allow the unit to rotate 180° and switch off when the adjustment screws are aligned again.
4. Unlock the 2 locking screws positioned on the drive boss using the extended wrench provided. The tilt angle can now be adjusted

using the short hex wrench. Release the adjusting screws and adjust until the desired angle is obtained and then lock off the adjusting screws. Lock off the locking screws on the drive boss.

Maintenance & Servicing

WARNING: Ensure the unit is disconnected from the mains electricity supply before attempting maintenance or servicing.

Periodically clean the instrument using a damp cloth and mild detergent solution. Do not use harsh or abrasive cleaning agents.

Any repairs or replacement of parts MUST be undertaken by suitably qualified personnel.

The following spares and accessories are available from your laboratory dealer.

Description	Cat Number
Spare securing bar (SI50)	RSI50
Blue Dyson non-slip mat (SI70 & 80)	KS0305
Spare platform tier (SI70 & 80)	SI700090

For a comprehensive list of parts required by service engineers conducting internal repairs and a service manual, please contact the Technical Service Department of Barloworld Scientific Ltd. quoting both the model and serial number.

Only spare parts supplied or specified by Barloworld Scientific Ltd. or its agents should be used. Fitting of non-approved parts may affect the performance and safety features designed into the instrument.

If in any doubt, please contact the Technical Service Department of Barloworld Scientific Ltd. or the point of sale.

Barloworld Scientific Ltd.
Stone, Staffordshire ST15 0SA
United Kingdom
Tel: +44 (0) 1785 812121
Fax: +44 (0) 1785 813748
e-mail equipment@barloworld-scientific.com
www.barloworld-scientific.com

5. Switch the motor back on and check operation.

Adjustment of Tilt Angle - SI80

1. Position the unit with its left side facing forward.

2. Remove the white plastic cap located in the centre at the bottom of the unit.

NOTE: The mains power needs to be connected during adjustment so care must be taken not to touch anything but the parts in the instruction below. All the live electrical parts are mounted on the front panel so are well away from the area of concern.

3. Turn the speed control to its slowest setting and, looking through the hole, watch for the black Allen screw to come into view. Stop the motor when it is directly facing you.
4. Using the tool provide with the rocker, unscrew the Allen screw by two turns. Remove the tool and rotate the motor by 180° and stop it when the silver headed slot screw is directly facing you.
5. Turn the screw clockwise until tight. This method increases tilt. Continue the above until the desired tilt is achieved.
6. To reduce the tilt, turn the screws in the opposite direction in the opposite sequence, i.e. undo the slot headed screw & tighten the Allen screw.

Warranty

Barloworld Scientific Ltd warrants this instrument to be free from defects in material and workmanship, when used under normal laboratory conditions, for a period of two (2) years. In the event of a justified claim Barloworld Scientific will replace any defective component or replace the unit free of charge.

This warranty does NOT apply if damage is caused by fire, accident, misuse, neglect, incorrect adjustment or repair, damage caused by incorrect installation, adaptation, modification, fitting of non approved parts or repair by unauthorised personnel.

Technical Specification (all models)

Temperature range	Ambient + 5° to 60°C
Temp. display resolution	0.1°C
Temperature accuracy	±0.5°C
Overall dimensions (wxdxh)	450 x 436 x 606mm
Internal height	321mm
Electrical supply	230V, 50Hz, 350W

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These products meet the relevant EC harmonised standards for radio frequency interference and may be expected not to interfere with, or be affected by, other equipment with similar qualifications. We cannot be sure that other equipment used in their vicinity will meet these standards

and we cannot guarantee that interference will not occur in practice. Where there is a possibility that injury, damage or loss might occur if equipment malfunctions due to radio frequency interference, or for general advice before use, please contact the Technical Service Department of Barloworld Scientific Ltd.

Declaration of Conformity

Model SI40, SI50, SI70, SI80 INCUBATORS

These products comply with the requirements of the EU Directives listed below:

89/336/EEC Electromagnetic Compatibility Directive, amended by 92/31/EEC and 93/68/EEC.

73/23/EEC Low Voltage Directive amended by 93/68/EEC

Compliance with the requirements of these Directives is claimed by meeting the following standards:

EN61326:1997, EN61000-3-2:1995, EN61000-3-2:1995, EN61000-4-4:1995, EN61000-4-5:1995, EN61000-4-6:1996, EN61000-4-11:1994, EN61010-1:1993 + Amendments A1:1992 and A2:1995, EN61010-2-010:1994

Compliance Certificates and Full Reports.

Ref: RETS0287/17/B, RETS0287/18/A

From "Epsilon Technical Services" an Independent Accredited Test House showing compliance to the above Standards are available on request

Signed:  (Mr C. R. Baggaley)

Date: 13/2/02

Authority: Technical Manager

for

Bibby Sterilin Ltd - Stone - Staffs - ST15 OSA - UK
Tel: +44 (0) 1785 812121 - Fax +44 (0) 1785 813748

INSPECTION REPORT

MODEL SI50, SI70, SI80

ELECTRICAL SAFETY

1. Earth continuity
2. Insulation
3. Flash test


FUNCTIONAL

1. Indicators
2. Speed / temperature
3. Visual acceptance

QUALITY CONTROL INSPECTOR

BZSj Scientific France Ltd

BP 79-77793 Nemours Cedex
France
Tel: +33 1 64 45 13 13
Fax: +33 1 64 45 13 00
e-mail: bibby@bibby-scientific.fr

BZSj Scientific Italia Ltd

Via Alcide de Gasperi 56
20077 Riozzo Di Cerro Al Lambro
Milano Italy
Tel: +392 983 8969
Fax: +392 982 30211
www.bibby-scientific.it

BZSj Scientific Ltd

Beacon Road Stone
Staffordshire ST15 0SA
United Kingdom
Tel: +44 (0)1785 812121
Fax: +44 (0)1785 813748
www.bibby-scientific.com

5j_RJRS 4 caž

350 Commerce Drive
Rochester NY 14623
USA
Tel: (800) 828-6595
Fax: (585) 334-0241
www.dynalabcorp.com

Afora S.A.

Calle Aribau 240
08006 Barcelona Spain
Tel: +343 93-306 98 00
Fax: +343 93-306 98 23
e-mail: marketing@afora.com
www.afora.com