



MANUFACTURED IN THE UK

Sterilin Ltd.

Angel Lane Aberbargoed Bargoed
Caerphilly CF81 9FW United Kingdom

Customer services:

Tel +44 (0)844 844 3737 Fax +44 (0)844 844 2373

Head office:

Tel +44 (0)1443 830 830 Fax +44 (0)1433 821 545

Registered No. 6512031

BS EN ISO 9001 registered company

TECHNICAL DATA SHEET

Product Group		POLYPROPYLENE DIPPAS®
Material Composition	Base	Polypropylene
	Caps	Polyethylene high and low density
	Handle	Blue Polypropylene, handle soldered to the container by ultrasound

Product Quality

Sterilisation By electron beam at 10Kgray (SAL ≈6)

Shelflife 3 years.

Product Details

Item Code	Description	Capacity (ml)	Colour	Handle Length (mm)	Materials Base/Cap	Sterility	Case Qty
195PPD	Dipper, screw cap	40	Blue	220	PP/PE	IRR	250
200PPD	Dipper, hinged cap	45	Blue	220	PP/PP	IRR	250
205PPD	Dipper, hinged cap	90	Blue	220	PP/PP	IRR	150
210PPD	Dipper, screw cap	125	Blue	220	PP/PE	IRR	100
215PPD	Dipper, screw cap	180	Blue	220	PP/PE	IRR	100



Quality in Life Sciences



IWAKI

