

SC6+ Colony counter

Improved version of SC6

Features:

- Energy saving LED lighting
- Multi voltage
- Connection to Printer
- Connection to PC
 - Software will be available to import to Excel
- Improved features at the same price!
- Available for shipment in North America June 2011




Printer
Connection

PC
Connection

SC6 Software 12	
Number of dishes.....	4
Dish 1.....	9
Dish 2.....	13
Dish 3.....	17
Dish 4.....	12
Mean count.....	12.8
Standard deviation n.....	2.9
Standard deviation n-1.....	3.3
CV n.....	+/- 22.4%
CV n-1.....	+/- 25.9%
Maximum count.....	17
Minimum count.....	9



- Replaced Fluorescent tube with LED's



Black Background is still available