



MANUFACTURED IN THE UK

Dynalab Corp. *Dynalon Labware*

350 Commerce Drive
Rochester, New York 14623 USA

Tel: 585-334-2060 800-828-6595

Fax: 585-334-0241

dynalabcorp.com

dynalon.com

TECHNICAL DATA SHEET

Product Group **MICROTITRATION PLATES**

Material Composition General Purpose Virgin Polystyrene
Colourant additive (black and white plates)

THE ABOVE MATERIALS ARE FDA APPROVED

Number of Wells 96

Well Capacity
U well 330ul
V well 310ul
Flat well 400ul

- Manufacture**
- 1 All Sterilin Microtitration Plates are manufactured in accordance with the Quality System Standard BS EN ISO9001:2000
 - 2 Verification of product clean room conditions in accordance with BS EN ISO 14644 'Environmental Cleanliness in enclosed spaces'
 - 3 Sterilin Polystyrene microtitration plates conform to the In-vitro diagnostic directive (98/79/EC)

Sterilisation

All Microtitration Plates sold as sterile are sterilised by gamma irradiation

Sterility is obtained through irradiation according to ISO 11137 (sterilisation of healthcare products- requirements for validation and routine control-radiation sterilisation) and EN552 (sterilisation of Medical Devices – validation and routine control of sterilisation by irradiation).

