

Levo Pipette Controller

Comfortable and easy to use



quality . value . service



Levo Pipette Controllers

- Comfortable – simple to use
- Precise pipetting – use with glass or plastic volumetric and graduated pipettes from 0.1 to 100mL
- Blow out feature
- Uses protective/replaceable syringe filter
- Easy to take apart and clean when needed
- 4 colors available for quick identification.



411125

Color
Length inches
Width inches
Price Each \$

-0001

blue
8
2
44.95

-0002

green
8
2
44.95

-0003

yellow
8
2
44.95

-0004

pink
8
2
44.95

Tel 800-334-7585

Fax 585-334-0241

www.dynalton.com

Dynalton
Labware
quality . value . service